




Joerns  
Oxford®

# Custom Sling Order Form

## Quickfit Glide

### STEP 1 Contact Details:

Company:	
Delivery Address:	
Tel:	
Email:	



Prior to requesting changes or if you're in any doubt, please consult this sling's standard specification on page 4.

### STEP 2 Core Sling Dimensions: Select a standard size or enter custom dimensions.

	<input type="checkbox"/> Small	<input type="checkbox"/> Medium	<input type="checkbox"/> Large	Custom
Dimension A	960mm	1030mm	1110mm	mm
Dimension B	660mm	640mm	650mm	mm
Dimension C	1270mm	1260mm	1290mm	mm

**Note:** Use page 3 to illustrate bespoke sizing. Dimension A is the mid-point of the sling body. All dimensions are subject to manufacturing tolerances.

### STEP 3 Material Options: Select the core material for your sling and include any bespoke requirements in the box provided.

<input type="checkbox"/> Polyester Solid	<input type="checkbox"/> Super-Soft
<input type="checkbox"/> Polyester Net	<input type="checkbox"/> 3D Spacer

Special Instructions:

**Note:** Some slings incorporate more than one fabric type - please consider this when requesting material changes. See page 4 for more information.

### STEP 4 Padding Options: Add padding to your sling for extra comfort. You can add a head support to your sling in STEP 6.

<input type="checkbox"/> 5mm HD Foam	<input type="checkbox"/> 5mm Soft Foam <input type="checkbox"/> 10mm Soft Foam	<input type="checkbox"/> 3mm Spacer	<input type="checkbox"/> Hollow Fibre
--------------------------------------	---	-------------------------------------	---------------------------------------

<input type="checkbox"/> Head Support
<input type="checkbox"/> Body
<input type="checkbox"/> Legs

Special Instructions:

### STEP 5 Inner Lining Options: A lining can be added to the inside of your sling for enhanced comfort.

<input type="checkbox"/> Micro Fleece	<input type="checkbox"/> Fixed (Stitched)	<input type="checkbox"/> Detachable (Velcro®)	<input type="checkbox"/> Head Support	<input type="checkbox"/> Body	<input type="checkbox"/> Legs
---------------------------------------	---	---	---------------------------------------	-------------------------------	-------------------------------

Special Instructions:

**STEP 6** **Head Support:** Use page 3 to help illustrate the shape and size of head support if required.☐ Add Head Support☐ Add Adjustment Straps☐ Add Stiffeners*Special Instructions:***STEP 7** **Straps, Loops, Clips & Handles:** Adjust strap lengths, number of loops, loop colours and positioning handles.

	Standard Strap Length	Custom Strap Length	Standard No. of Loops	No. of Extra Loops?
Shoulder	400mm	mm	4	
Hip	N/A	mm	N/A	
Leg	400mm	mm	3	

*Special Instructions: (e.g. quantity and placement of positioning handles)*☐ Add Positioning Handle(s)☐ Omit Positioning Handle(s)☐ Omit Modesty Loop☐ Clip Attachments (Securi3®)*When choosing clips, please specify the length of strap required from the sling body to the clip (required for both shoulder and leg attachment points).*

**Loop Colour Options:** Black, blue, green, grey, orange, red, white and yellow. Please confirm selection/sequence in the special instructions box above. See page 3 for dimensional illustrations for standard strap lengths and loop positions.

**STEP 8** **Body Belt:** Add or alter a body belt to enhance support and security during transfer. Use page 3 to help illustrate your needs.☐ Buckle**Belt Length****Belt Width***Special Instructions: (e.g. padding, elastic, etc)*☐ Velcro®

mm

mm

**STEP 9** **Other Options:** Use page 3 to help illustrate your bespoke requirements.☐ Add Positioning Pocket☐ Add Centre Markers*Special Instructions:***STEP 10** **Additional Requirements****STEP 11** **Confirmation**☐ I am happy with and hereby give permission to proceed with the custom sling order defined on this form.☐ Please call me to discuss/confirm this order and/or price prior to manufacture.**Phone No:****Customer PO Number:****Order Quantity:****Sales Order Number:****Ordered By:****Signed:****Date:**

**Please Note:** Joerns Healthcare reserve the right to refuse a special order if deemed unsafe. In such event, we will contact you to discuss the matter further.

**Please return your completed form to:**


**orders@joerns.co.uk**

**Please attach any supporting illustrations when submitting your form.**

**(T)** +44 (0)344 811 1158 **(W)** www.joerns.co.uk

# Illustrative Guidance

Please use this page to help illustrate or clarify any bespoke requirements



# Standard Specification Breakdown

Please consider the following information when requesting changes to your standard sling.

	<b>Sling Model</b>
	<b>Oxford Quickfit Glide</b>
<b>Material</b>	
Core Material	Polyester Solid (Grey)
Spine Overlay	Rip-stop Nylon (Black)
Leg Overlay	Rip-stop Nylon (Black)
<b>Padding</b>	
Head Support	N/A
Body	5mm Soft Foam to Spine
Leg Section	5mm HD Foam
<b>Head Support</b>	
Included	N/A
Integrated Stiffeners	N/A
Adjustment Straps	N/A
<b>Loops, Clips &amp; Handles</b>	
Sling Attachment Type	Loops
Positioning Handles	2
Modesty Loop	✓
<b>Body Belt</b>	
Included	N/A
Attachment Type	N/A
<b>Other Options</b>	
Centre Markers	N/A
Positioning Pocket	N/A

You do not need to submit this page, it is for your reference only.